



APEX
MOBILITY

THOMASHILFEN BROCHURE

apexmobility.com.au | 1300 212 192



The APEX Mobility Product Range

The APEX Mobility team continually research new products and global innovation, which is why we can provide the very best assistive technology equipment available in the Australian market today. Before any products are added to the APEX Mobility range, they are tested and trialed to make sure they meet all functional and quality requirements, plus regulatory standards for safety.

Products featured in this brochure are part of our dynamic range of quality assistive technology equipment. Most of these products are scripted/prescribed, which means they are customised to suit the specific needs of individual clients.



Thomashilfen

“

*Our goal is your health, well-being
and positive development.*

”

Pushchairs

Jazz Advantage S
Swiftly

3 – 5

Car Seats

Commander
Defender

6 – 9



Jazz Advantage S

Sitting for long periods can be tiring, but sometimes it cannot be avoided. The Jazz Advantage S rehab stroller makes life easier with its many setting options, adaptability, and super soft upholstery for cosy comfort when sitting for extended periods.

Key Features:

- Soft suspension absorbs vibrations for a more enjoyable ride without compromising steering.
- Easy one hand adjustment to tilt the seat for a rest position or active sitting.
- Forward and rear facing options so the child can look in the direction of travel or to the carer.
- No tool adjustment features for individual and therapeutic adaption including the seat width and dept, back height and footrest.
- 130° recline of seat for lying position if needed.
- Handy storage space including a washable basket. A large tray (optional accessory) can be attached as a vent tray for medical devices.
- Thick cushions which are covered with breathable and washable fabric.
- Choice of 5 upholstery styles/colours.

Specifications:

Jazz Advantage S	Size 1	Size 2
Seat Depth	19-30cm (7.5"-11.8")	27-38cm (10.6"-15")
Seat Width	19-31cm (7.5"-12.2")	25-37cm (9.8"-14.6")
Seat Height	47cm (18.5")	46cm (18.1")
Back Height	41.5-68cm (16.3"-26.7")	55-74cm (21.7"-29.1")
Hip Angle	90° to 130°	90° to 130°
Seat Tilt-In-Space	10° to +35°	10° to +35°
Overall Dimensions (l x w x h)	99.5 x 59.5 x 114cm 39.2 x 23.4 x 44.9"	107 x 69 x 120cm 42.1 x 27.2 x 47.2"
Weight (seat unit)	9kg (19.8 lbs)	10.4kg (22.9 lbs)
Weight (chassis)	9.8kg (21.6 lbs)	10.9kg (24 lbs)
Max. user weight	35kg (77.2 lbs)	40kg (88.2 lbs)



Accessories:



Vent Tray



Grip Rail



Abduction Block



Lateral Trunk Supports



Hi-Lo Indoor Base



Accessory Bag



Sun Canopy



Tray

Upholstery Colours:



Swiftly

Swiftly is a lightweight, stable and comfortable paediatric pushchair with a fixed seat. Swiftly has a simple folding system for fast transportation and can withstand the many stresses of day-to-day life. Available in two sizes for growing children.

Key Features:

- Back angle can be adjusted for positioning and relaxed seating.
- Seat, back and footrest provide safe and firm support.
- Seat depth adjustable without tools, even when the child is sitting in the chair.
- H-Belt harness with padded harness as standard feature.
- Seat minimiser for smaller children. Grows with the child.
- Size 2 Swiftly designed for children weighing up to 50kg (110 lbs).
- Foldable and height adjustable footrest, which may assist children to get in and out of the stroller independently.



Specifications:

Swiftly	Size 1	Size 2
Seat Depth	22 (17cm) – 28.5 cm	37.5cm (28cm)
Seat Width	34cm (23cm)	39.5cm (29.5cm)
Back Height	62cm / 24.4"	70cm / 27.6"
Overall Dimensions (l x w x h)	61 x 102 x 98cm 24 x 40.2 x 38.6"	67 x 113 x 98cm 26.4 x 44.5 x 38.6"
Product weight	12.4kg (27.3 lbs)	18.2kg (40.1 lbs)
Max. user weight	35kg (77.2 lbs)	50kg (110.2 lbs)

Accessories:



Lateral Trunk Support



Grip Rail



Abductor Block/Pommel



Headrest




Tip Assist



Sun Canopy



Basket



Tray



Commander

Commander is a highly adaptive positioning car seat for children with special needs. Especially for children with little upper body control, the fixed lateral trunk supports offer this extra support and protection – the supports have proven themselves over many years of use in rehab strollers. A huge selection of accessories and easy handling make the adjustment perfect.

Key Features:

- The unique, modular combination of swivel base and tilting base can be supplemented with a variety of accessories as required.
- With the integrated seat tilt-in-space +3° a resting position is quickly adjusted.
- The breathable cover can be completely removed and washed.
- Release lever (draw latch) for easy handling and easy to reach. It can be mounted either on the left or right front of the swivel base.
- Comfortable handling - 5-point harness automatically adjusts when growing with the child.
- Integrated 9-way height-adjustable headrest including additional padding.
- Optional use of modular function plates with rotation and/or tilt function and further accessories.
- Handy cup holder (right or left applicable).

Specifications:

Commander	Measurement
Seat depth	34 cm / 13.4"
Seat width, front (with supports)	33.5 cm (22 cm) / 13.2" (8.7")
Seat width, rear (with supports)	28.5 cm (17 cm) / 11.2" (6.7")
Back height (seat)	55 - 75 cm / 21.7 - 29.5"
Shoulder height (belt guide)	25 - 44 cm / 9.8 - 17.3"
Seat tilt-in-space (basic model)	3°
Seat tilt-in-space (with tilting plate)	5.5° / 11° / 16°
Lower leg length (Footrest directly attached to adapter plate)	23 - 40 cm / 9.1 - 15.7"
Lower leg length (Footrest attached to swivel base + tilting base)	28 - 45 cm / 11 - 17.7"
Overall dimensions (w x h x d)	55 x 65 x 45 cm / 21.7 x 25.6 x 17.7"
Weight of seat	approx. 8.5 kg / 18.7 lb
Max. user weight	15 - 36 kg / 33.1 - 79.4 lb

Accessories:



Abduction block and hip supports support a stable sitting position.



Lateral trunk supports (metal core) can be positioned precisely on the back plate with hook-and-loop-fasteners for individual adjustment.



Crotch pad (large) protects the sensitive belly of the child, **SEAT WEDGE (BELOW)** allows +10° seat tilt.



Seat wedge (inside) supports an upright pelvic position, **LATERAL TRUNK SUPPORTS (METAL CORE)** for optimum support of the upper body.



Table promotes the straightening up and stable posture of the child, provides the child with all-round support and safety.



Anti-escape vest prevents the child from wriggling out of the belt.



Car seat protection mat with back protects the seat upholstery and offers 3 practical compartments for toys at the front.

Upholstery Colours:





Defender

Defender is a 2-in-1 car seat solution that helps families enjoy safe travelling. A car seat & booster, the Defender offers children with special needs 360° of impact protection. The head, trunk and thighs supports, combined with unique accessories offer many positioning options. The seat allows for growth while remaining narrow enough for a second car seat as well as another child to fit in the back seat.

Key Features:

- Head, trunk and thighs positioning supports inclusive.
- 9-position height adjustable headrest.
- Lateral side supports are removable to adjust seat width for plenty of room for the child to grow with the seat.
- No-threaded 5-point harness provides additional support in Booster Mode.
- Ultra plush seat pad made with durable premium fabrics and extra foam for maximum comfort. It is removable for machine washing.
- Cup/snack holder, chest clip and soft harness covers.
- Optional swivel base and footrest.

Specifications:

Defender	Measurement
Seat depth	34.3 cm / 13.5"
Seat width*	21.6 – 30.5 cm / 8.5 – 12"
Seat to head**	60 – 77.5 cm / 23.5 – 30.5"
Shoulder height (child)***	25.4 – 49.5 cm / 10 – 19.5"
Lower leg length	(18 cm) 28 – 43 cm / (7") 11 – 17"
Overall dimensions (w x h x d)	46.3 x 67.3 x 50.8 cm / 18.25 x 26.5 x 20"
Weight of seat	10 kg / 22 lb
Harness Mode	Measurement
Weight of child	10 – 29.5 kg / 22 – 65 lb
Height of child	68.6 – 144.8 cm / 27 – 57"
Booster Mode	Measurement
Weight of child	13.6 – 50 kg / 30 – 110 lb
Height of child	86.4 – 144.8 cm / 34 – 57"

* To reach max. seat width remove additional hip guide pads from additional seat pad.

** To reach max. back height remove additional seat pad and foam pad under seat cushion.

*** To reach max. shoulder height remove additional seat pad and foam pad under seat cushion.

Accessories:



Footrest, short with footrest adapter and footrest pad



Footrest, long with footrest adapter



Seat wedge, below (+ 10°)



Table = playtray and additional impact shield



Abduction block



Crotch pad, large

Product Demonstrations & Trials

APEX Mobility does not have demo units available. In many cases, Dealers are offered subsidised pricing for units that they will hold for showroom display or trial use only.

Dealers of Assistive Technology

Contact your APEX Mobility account manager to discuss your needs.

Therapists

To arrange a trial of any of this range, contact the following dealer:

Authorised Dealer



 1300 212 192  sales@apexmobility.com.au

 apexmobility.com.au

V2 September 2022